



MINNTECH[®]
renal systems

A Minntech Corporation Business Group

Product Guide


Reprocessing
and Dialysis
Supplies

U.S. ORDERS: 1-800-328-3340
www.minntech.com/renal

VOLUME 6 | ISSUE 2

Reprocessing Accessories

Dialyzer Caps and Connectors



P/N	Description	Package
78196-010	Dialysate Port Caps (Red) - A	100/bag
78196-050	Dialysate Port Caps (Red) - A	500/bag
78398-367	Ventable Dialysate Port Caps, 1.5" - B	50/bag
78208-000	Ventable Dialysate Port Caps, 2.25" - B	50/bag
78197-010	Blood Port Caps - C	100/bag
78197-100	Blood Port Caps - C	1,000/bag
78397-704	ISO Fitting Blood Port Caps - D	100/bag
78398-366	ISO Fitting Blood Port Caps - D	1,000/bag
78397-713	Reprocessing Connectors (Barrier Adapter) - E	10/bag
78397-699	Reprocessing Connectors (Barrier Adapter) - E	100/bag
78398-046	Reprocessing Connector Caps - F	100/bag



Caps Container

RENALIN Cold Sterilant resistant, graduated, 2 quart plastic container with lid. Caps Container can hold dialysate port caps, reprocessing connectors or blood port caps for easy access.

P/N	Description	Package
47049-142	Caps Container, Square	Each

Reprocessing Supplies to Meet Your Needs

Renalog[®] TM –



A Comprehensive Training Program for Dialyzer Reprocessing

Call Clinical Services at 1-800-328-3345 for details.

Terms: F.O.B. Minneapolis, MN. Effective June 1, 2011. Call your Minntech Representative for current pricing information.

20 | For U.S. orders call 1-800-328-3340



Reprocessing Labels & Accessories



Label Kit, Renalog RM
Supplies to support the Zebra® TLP 2844 label printer, includes 4 label rolls, 4 dye transfer ribbons, and a printer cleaning card.

P/N	Description	Package
78398-179	Label Kit, RENALOG RM	Each



Labels, Computer Form-Fed
For use with RENALOG Software using dot matrix printers. Individually perforated. 5,000 per pack.

P/N	Description	Package
40080-025	Form Fed Labels for RENALOG	Pack

Labels, Dialyzer



P/N	Description	Package
40080-030	Dialyzer Barcode Label Protector	100/roll
40080-035	Dialyzer Barcode Label Protector (Similar Name)	100/roll
67655-001	DRS-002 White (Opaque) Labels	100/roll



Labels, Hazard I.D.
Identify workplace bottles and containers with the hazard information workers need to remain safe. This economical label lets you write in the chemical name, manufacturer, date and hazard rating. Perforated on a roll of 720 labels.

P/N	Description	Package
78398-058	Hazard I.D. Labels	Roll



Labels, Similar Name
These RENALIN Cold Sterilant resistant, bright orange labels are ideal for identifying those patients with similar names. Packaged in a convenient dispenser box. 120 per roll.

P/N	Description	Package
78398-039	Similar Name Labels	Roll



Labels, RENALIN 1%
Used when making a 1% RENALIN solution, to identify the person who filled the container and the date it was filled. 100 per roll.

P/N	Description	Package
67198-251	RENALIN 1% Label	Roll

Questions? ...we have answers!

- Free patient education material
- Resource CD
- Renew - dialyzer reprocessing newsletter

Visit www.minntech.com



Terms: F.O.B. Minneapolis, MN. Effective June 1, 2011. Call your Minntech Representative for current pricing information.

Test Kits



Test Strips, PERASSAY® 500 Peracetic Acid **
A dip and read, go/no-go measurement of peracetic acid at a 500 ppm level to confirm the presence of RENALIN solution in reprocessed dialyzers after storage and prior to rinsing for clinical use. 100 per vial.

P/N	Description	Package
78378-000	PERASSAY 500 Peracetic Acid Test Strips	Vial

** Indicates components of the RENALIN PR Test System™



Test Strips, GUARDIAN™ SERIM® Chlorine Residual
Tests for residual levels of chlorine in rinse water. Tests can be completed in 5 seconds for a qualitative measure and 30 seconds for a semi-quantitative measure. Measures down to 0.5 ppm. Six bottles of 100 test strips per box.

P/N	Description	Package
78400-582	SERIM Residual Chlorine Test Strips	Box



Test Strips, RENALIN Residual **
Tests for a quantitative measure of RENALIN concentration after rinse-out. This test is sensitive to one part per million. 100 per vial.

P/N	Description	Package
78198-000	RENALIN Residual Test Strips	Vial

** Indicates components of the RENALIN PR Test System™



Test Strips, GUARDIAN™ SERIM® Hardness
Takes a semi-quantitative measurement of total hardness in the water. Measures 0-10 ppm. Six bottles of 50 test strips per box.

P/N	Description	Package
78400-581	SERIM Hardness Test Strips	Box



Test Strips, ACTRIL Indicator
The ACTRIL Indicator Test Strip indicates the presence of ACTRIL Cold Sterilant within minutes. 500 per vial.

P/N	Description	Package
78259-000	ACTRIL Indicator Test Strips	Vial



RENATEST Deluxe Pump Kit RH 1100
The RENATEST Vapor Detection System aids in maintaining staff safety and compliance with OSHA and CMS regulations.

P/N	Description	Package
67757-040	RENATEST Deluxe Pump Kit RH 1100	Each



Test Strips, ACTRIL® Residual
The ACTRIL Residual Test Strip confirms that ACTRIL Cold Sterilant concentrations have been reduced to safe levels during rinse procedures. 100 per vial.

P/N	Description	Package
78258-000	ACTRIL Residual Test Strips	Vial



Detector Tube Supplies for the RENATEST Kit
10 per case.

P/N	Description	Package
67757-031	RH-1010 Hydrogen Peroxide Detector Tubes (0.1-3 ppm)	Case
67757-032	RH-1020 Acetic Acid Detector Tubes (5-80 ppm)	Case



Deluxe Tube Cutter
RENATEST, replacement Deluxe Tube Cutter.

P/N	Description	Package
78398-515	Deluxe Tube Cutter	Each

*Visit the Minntech Website
www.minntech.com*

Terms: F.O.B. Minneapolis, MN. Effective June 1, 2011. Call your Minntech Representative for current pricing information.

22 | For U.S. orders call 1-800-328-3340





Test Strips, SteriChek® Bicarb pH Reagent

These strips provide a convenient and accurate means of measuring the pH of bicarbonate concentrates and solutions of bicarbonate and acid. The test strips are supplied ready-to-use. The indicator pad changes color relative to the level of hydrogen ion concentration in bicarbonate buffer when the strip is used according to the directions. Kit of 5 bottles of 100 strips.

P/N	Description	Package
78399-261	SteriChek Bicarb pH Reagent Strips	Box



Low-Range Water Hardness Test Strips

SteriChek® Sensitive Low-Range Water Hardness Strips have color blocks at levels of 0, 5, 10, and 20 ppm (mg/l) or 0, 0.3, 0.6, and 1.2 ppg. The test strips are supplied ready-to-use. The indicator pad changes color relative to the level of calcium and magnesium present when the strip is used according to the directions. 50 strips per bottle. 6 bottles per box.

P/N	Description	Package
78398-595	Water Hardness Test Strips	Box



Test Strips, Serim® HiSense Ultra 0.1 Chlorine

A quick, convenient and accurate means for testing low levels of Total Chlorine (total chloramine) in feed water and Residual Chlorine in rinse water. Measures down to 0.1 ppm (mg/l) Chlorine. Meets AAMI standards. FDA cleared. Easy-to-distinguish color blocks. Five 100 strip vials per box.

P/N	Description	Package
78400-290	Serim HiSense Chlorine Test Strips	Box



Test Strips, pH 0 - 14

Universal pH test strips that test the full 0-14 pH range. The pH indicator strips are for the rapid measurement of the pH of a water sample. They can also be used for monitoring the pH of a sample during analysis to avoid errors due to weak buffered solutions. Indicator dyes are chemically bonded to the strip to avoid the possibility of color bleeding in solutions. 100 per box.

P/N	Description	Package
78398-045	pH Test Strips 0-14	Box



Total Chlorine DPD Test Kit

Measures low levels of total chlorine (i.e. total chloramines plus free chlorine in feed water used to prepare dialysate).

The test kit will provide quantitative results in increments of 0.1 ppm utilizing the recognized standard method of analysis of the Association for the Advancement of Medical Instrumentation (AAMI). The kit includes reagent, apparatus, and quality control measure for monitoring results. 100 DPD powder pillows, two sample tubes, and one color comparator wheel per kit.

P/N	Description	Package
78398-596	Chlorine DPD Test Kit	Kit



WaterCheck™ CP

Easy to use test strip can be used as a residual or potency test for chlorine. Detects 0 to 5000 ppm chlorine. Easy to read color blocks at 0, .5, 1ppm chlorine, 1% bleach and 10% bleach.

P/N	Description	Package
78400-296	WaterCheck™ CP Bleach Test Strips	1 Vial/100



Total Chlorine DPD Refill Kit

Refill kit for the Total Chlorine DPD Test Kit. Kit includes 500 powder pillows per box.

P/N	Description	Package
78398-597	Chlorine DPD Refill Kit	Box

A Variety of Strips Available to Meet Your Testing Needs

Terms: F.O.B. Minneapolis, MN. Effective June 1, 2011. Call your Minntech Representative for current pricing information.

renalin® 100 Cold Sterilant

The gold standard in dialyzer reprocessing

- Effective
- Environmentally friendly
- No pre-dilution
- Vented caps prevent pressure build-up
- Non carcinogenic
- Ready-to-use



RENALIN 100 Cold Sterilant Concentrate

For dialyzer reprocessing with the Renatron® 100 Series system. Each 1 gallon container holds three liters of RENALIN 100 concentrate. To prevent excessive pressure build-up, all concentrate containers are provided with vented caps. These caps should not be altered or replaced. All boxes must be stored in an upright position. Four 1 gallon containers per case.

P/N	Description	Package
78397-844	RENALIN 100 Cold Sterilant Concentrate	Case
78397-845	RENALIN 100 Cold Sterilant Concentrate with 100 RENALIN Residual & 100 Perassay® 500 Test Strips	Case
78397-846	RENALIN 100 Cold Sterilant Concentrate with 300 RENALIN Residual & 300 Perassay® 500 Test Strips	Case

Terms: F.O.B. Minneapolis, MN. Effective June 1, 2011. Call your Minntech Representative for current pricing information.

Sterilants, Disinfectants, & Cleaners



ACTRIL® Cold Sterilant

ACTRIL Cold Sterilant provides proven and effective surface disinfection while safeguarding your patients and staff. Four 1 gallon containers per case.

P/N	Description	Package
78270-000	ACTRIL Cold Sterilant System with 100 Residual & 100 Indicator Test Strips	Case
78337-000	ACTRIL Cold Sterilant	Case



Formula 409® - For Institutional Use

The only Formula 409 formulation validated as a system cleaner agent for the RENATRON system.

P/N	Description	Package
93250-092	Formula 409	1 gallon



Machine Tag, RENALIN

These 4" x 3.5" Florescent Orange cards attach conveniently to the Hansen connector of the RENATRON or delivery system.

P/N	Description	Package
67198-444	Machine Tag, RENALIN	Each



Machine Tag, Formula 409®

These 4" x 3.5" Florescent Pink cards attach conveniently to the Hansen connector of the RENATRON or delivery system.

P/N	Description	Package
67198-445	Machine Tag, Formula 409	Each



CleanCide® Ready-To-Use

Ready to use formula meets OSHA requirements for bloodborne pathogens (EPA Registered). "One-step" cleaning and disinfecting process is effective against: MRSA, E. Coli, VRE, Hepatitis B and more. Does not contain phenol, quats, alcohol or bleach. Active Ingredient: Citric Acid. 12 quarts per case.

P/N	Description	Package
78400-271	Cleancide RTU	Case



CleanCide™ Germicidal Wipes

The wipes provide a "one-step" cleaning and disinfecting process effective against: MRSA, Escherichia Coli, VRE, Hepatitis B and more. CleanCide wipes meet OSHA requirements for bloodborne pathogens and are EPA Registered. CleanCide does not contain phenol, quats, alcohol or bleach. The Active Ingredient is Citric Acid. Each case contains 12 containers.

P/N	Description	Package
78400-272	CLEANCIDE Wipes - 6" x 6"	Case
78400-273	CLEANCIDE Wipes - 10" x 12"	Case



Chemflex®

Refillable closed container with dry wipes validated for use with 1% and 10% bleach, CleanCide® and Actril®. Eliminates stability issues with bleach and other disinfectants. Reduces the risk of cross-contamination, spills and vapor exposure. Low lint and Cellulose free. Six rolls per case.

P/N	Description	Package
78400-287	CHEMFLEX Wipes	Case



DrainGel™

Professional strength bacterial treatment that will safely and quickly biodegrade biofilm. Its powerful formula will also eliminate drain odors and small-fly infestations. Application is quick and easy. DrainGel can be poured down the sink and floor drains or applied using the DrainGel Sprayer. Safe for all types of plumbing.

P/N	Description	Package
78398-047	DrainGel	4 gal/case
78398-048	DrainGel Sprayer	Each

Terms: F.O.B. Minneapolis, MN. Effective June 1, 2011. Call your Minntech Representative for current pricing information.

On Demand® Bleach Autodilution Systems




ON DEMAND®
Autodilution System
for 1% Bleach



ON DEMAND®
Autodilution System
for 10% Bleach



WaterCheck™ CP
Bleach Test Strips



ON DEMAND® Autodilution System for 1% and 10% Bleach

The ON DEMAND AUTODILUTION SYSTEM was developed in response to a need for a safe, accurate and labor-saving alternative to manual preparation of disinfectant solution. At the touch of a finger, the unique venturi system proportions and dispenses a consistent, continuous supply of 1% or 10% bleach solution. This easy-to-install device requires no electrical connection or special electronic systems and can be used in minutes. Validation studies of the patented design of the ON DEMAND AUTODILUTION SYSTEM demonstrated consistent and accurate performance in the preparation of 1% and 10% bleach solution.

Verify solution Strength with WaterCheck™ CP Bleach Test Strips

P/N	Description	Package
78400-240	Bleach 1% Autodilution System	Each
78400-241	Bleach 10% Autodilution System	Each
78399-901	Water Hose	Each
78400-296	WaterCheck™ CP Bleach Test Strips	1 Vial/100

Terms: F.O.B. Minneapolis, MN. Effective June 1, 2011. Call your Minntech Representative for current pricing information.