



LONG TERM HAEMODIALYSIS CATHETER

Kflow® Epic® RETRO HAEMO

LONG TERM HAEMO Kflow® Epic® RETRO HAEMO


LONG TERM

ORDERING INFORMATION

Code	Size	Length* Tip to cuff	Length* Total	Description	Priming Vol A	Priming Vol V	Box Qty
KFE-RHDL-1216-K	12	16	21	Long term retro haemo stepped tip dual lumen catheter kit	0.90	1.00	10
KFE-RHDL-1219-K	12	19	24	Long term retro haemo stepped tip dual lumen catheter kit	1.10	1.20	10
KFE-RHDL-1419-K	14	19	24	Long term retro haemo stepped tip dual lumen catheter kit	1.35	1.45	10
KFE-RHDL-1423-K	14	23	28	Long term retro haemo stepped tip dual lumen catheter kit	1.50	1.60	10
KFE-RHDL-1427-K	14	27	32	Long term retro haemo stepped tip dual lumen catheter kit	1.60	1.65	10
KFE-RHDL-1431-K	14	31	36	Long term retro haemo stepped tip dual lumen catheter kit	1.70	1.75	10
KFE-RHDL-1435-K	14	35	40	Long term retro haemo stepped tip dual lumen catheter kit	1.80	1.90	10
KFE-RHDL-1450-K	14	50	55	Long term retro haemo stepped tip dual lumen catheter kit	2.55	2.60	10
KFE-RHDL-1619-K	16	19	24	Long term retro haemo stepped tip dual lumen catheter kit	1.50	1.60	10
KFE-RHDL-1623-K	16	23	28	Long term retro haemo stepped tip dual lumen catheter kit	1.70	1.80	10
KFE-RHDL-1627-K	16	27	32	Long term retro haemo stepped tip dual lumen catheter kit	1.90	2.00	10
KFE-RHDL-1631-K	16	31	36	Long term retro haemo stepped tip dual lumen catheter kit	2.00	2.10	10
KFE-RHDL-1635-K	16	35	40	Long term retro haemo stepped tip dual lumen catheter kit	2.45	2.60	10
KFE-RHDL-1650-K	16	50	55	Long term retro haemo stepped tip dual lumen catheter kit	3.10	3.20	10
KFE-RHDL-1216-PCK	12	16	21	Long term retro haemo stepped tip dual lumen pre curved catheter kit	0.90	1.00	10
KFE-RHDL-1219-PCK	12	19	24	Long term retro haemo stepped tip dual lumen pre curved catheter kit	1.10	1.20	10
KFE-RHDL-1419-PCK	14	19	24	Long term retro haemo stepped tip dual lumen pre curved catheter kit	1.35	1.45	10
KFE-RHDL-1423-PCK	14	23	28	Long term retro haemo stepped tip dual lumen pre curved catheter kit	1.50	1.60	10
KFE-RHDL-1427-PCK	14	27	32	Long term retro haemo stepped tip dual lumen pre curved catheter kit	1.60	1.65	10
KFE-RHDL-1431-PCK	14	31	36	Long term retro haemo stepped tip dual lumen pre curved catheter kit	1.70	1.75	10
KFE-RHDL-1435-PCK	14	35	40	Long term retro haemo stepped tip dual lumen pre curved catheter kit	1.80	1.90	10
KFE-RHDL-1450-PCK	14	50	55	Long term retro haemo stepped tip dual lumen pre curved catheter kit	2.55	2.60	10
KFE-RHDL-1619-PCK	16	19	24	Long term retro haemo stepped tip dual lumen pre curved catheter kit	1.50	1.60	10
KFE-RHDL-1623-PCK	16	23	28	Long term retro haemo stepped tip dual lumen pre curved catheter kit	1.70	1.80	10
KFE-RHDL-1627-PCK	16	27	32	Long term retro haemo stepped tip dual lumen pre curved catheter kit	1.90	2.00	10
KFE-RHDL-1631-PCK	16	31	36	Long term retro haemo stepped tip dual lumen pre curved catheter kit	2.00	2.10	10
KFE-RHDL-1635-PCK	16	35	40	Long term retro haemo stepped tip dual lumen pre curved catheter kit	2.45	2.60	10
KFE-RHDL-1650-PCK	16	50	55	Long term retro haemo stepped tip dual lumen pre curved catheter kit	3.10	3.20	10

*Length in cm

Tests done on extracorporeal circuit with fluid having same viscosity (3.5 cP at 37°C) of blood.



KIT CONTENTS

1 x Catheter	2 x Vessel dilator	2 x Injection caps
1 x 18g x 7cm introducer needle with side valve	1 x Peel away sheath introducer	1 x Sutureless securement device
1 x 0.035" x 70cm J tip nitinol guidewire	1 x Tunnelling stylet	2 x Wound dressing
	1 x 5ml syringe	1 x Patient identification sticker
	1 x Scalpel # 11	

Please contact your local representative for more details.

©2012 Kimal plc. All rights reserved. No part of this document may be reproduced in any form or by any means without permission in writing from Kimal plc. Specification: Kimal plc reserves the right to change the specification of its products at any time. Registered address: Kimal plc, Arundel Road, Uxbridge, UB8 2SA.

KIMAL
Making a Difference

Kimal UK
Unit 401, Pointon Way, Stonebridge
Cross Business Park, Droitwich Spa,
Worcestershire WR9 0LW
Tel: +44 (0)845 437 9542
Fax: +44 (0)845 437 9544
www.kimal.co.uk

Kimal Deutschland GmbH
Benediktusstrasse 88
40549 Düsseldorf
Tel: +49(0)211-58588125
Fax: +49(0)211-58588126
www.kimal.co.uk

Kimal International Healthcare FZE
SAIF - Zone
PO Box 8532, Sharjah
United Arab Emirates
Tel: +00 971 67451597
Fax: +00 971 65588146
mea@kimal.co.uk

Rev No.2

L-132B-0003-A1

ORDERING INFORMATION

Code	Description	Box Qty
UNIFIX	Unifix Adhesive sutureless securement device for catheter fixation and anchoring	50
KFE-RPK-AV	Repair extension assembly each box contains x3 A (red) & x3 V (blue) extension assemblies	6



KFE-RPK-AV



Unifix

Adhesive sutureless securement device

The Kflow Epic Retro Haemo catheter range has been developed for long term haemodialysis access. The biocompatible carbothane lumen construction is key to providing excellent tensile strength, flexibility and superior kink resistance. Available in straight and pre curved configurations in a variety of French sizes and lengths to best meet the needs of patients and clinicians.

FEATURES & BENEFITS

- Available in a variety of fr sizes to enable clinicians to choose the optimum size for the patient.
- Specifically designed tip configuration to reduce recirculation.
- Round catheter permits easy insertion and supports traditional modified Seldinger insertion technique
- Catheters available as straight or pre curved configuration. A pre curved catheter often used for jugular approach to make catheter more comfortable for the patient.
- Catheters suitable for standard or retrograde tunnelling. Female luer and extension line assemblies are provided to be assembled to catheter extensions after retrograde tunnelling. These can also be used for repairing catheter extensions.
- Nitinol guidewire as standard provides excellent kink resistance, and comes pre loaded in holder with thumb feeder capability to help clinician insert guidewire smoothly into the vessel.
- Round lumens enable blood flow rates of 350-400 mL/min while maintaining low arterial and venous pressures.
- Radiopaque carbothane material is solvent resistant withstanding all major iodine and alcohol based cleaning solutions.
- Priming Volume printed on catheter extensions to assist clinicians with catheter maintenance.
- Colour coded clamps and luers ensures easy lumen identification.



A Easily identifiable colour coded clamps and printed extensions



B Robust colour coded luers



C Cuff
To promote optimum tissue granulation and catheter securement



D Soft stepped catheter tip design aids smooth insertion and reduces recirculation rates



Extension assembly



Flexible kink resistant Nitinol guidewire

Unifix product here



Unifix
Adhesive sutureless securement device